NSK

Multi-Control Panel for built-in system

Multi Pad

Multi Pad for NLX / iMD

OPERATION MANUAL

Please read this Operation Manual carefully before use, and file for future reference.

Thank you for purchasing Multi Pad.

Read this Operation Manual carefully before use for operation instructions and care and maintenance guidelines.

Keep this Operation Manual for future reference.

Contents

⚠ CAUTIONS FOR HANDLING AND OPERATION	2
1. FEATURES	6
2. SPECIFICATIONS	6
3. NAME OF EACH PART	6
4. START UP	6
5. DESCRIPTION FOR DISPLAY	7
5-1. MOTOR	7
5-2. SCALER	8
5-3. SELECT THE INSTRUMENT	8
6. SETTINGS	9
6-1. MOTOR	9
6-2. SCALER	9
7. USEFUL FUNCTIONS	9
7-1. MOTOR	9
8. ERROR CODES	10
8-1. NLX nano MOTOR	10
8-2. NBX MOTOR	11
8-3. SCALER	11
9. OTHER SETTINGS	12
9-1. MOTOR SETTINGS	12
9-2. SCALER SETTINGS	15
10. DISPOSING PRODUCT	15
11. WARRANTY	15

CAUTIONS FOR HANDLING AND OPERATION

CONNECT ONLY NLX(nano) / NBX MOTOR. DO NOT CONNECT OTHER MOTOR.

- Read these safety cautions thoroughly before use and operate the product properly.
- These indicators will show you how to operate the product safely and prevent danger to you or others. They are classified by degree and/or severity of danger. All contents relating to safety should be observed.

Classification	Degree and severity of danger or damage
 DANGER	Provides an instruction where death or serious injury may occur.
 WARNING	Provides an instruction where personal injury or physical damage may occur.
 CAUTION	Provides an instruction where minor to medium injury or physical damage may occur.
N O T I C E	Provides an instruction that should be observed for safety reasons.

DANGER

- Do not attempt to disassemble the product or tamper with the mechanism. It may cause an electric shock or fire.

WARNING

- Do not use for implant surgery.
- The product is designed only for clinical dental use by qualified personnel.

CAUTION

- The Multi Pad operation panel is not water proof and should not be exposed to water.
- Do not use the Multi Pad as a handle when you want to move the dental unit.
- Do not pull the motor cord with excessive force.
- Be sure to press the buttons with a finger.
- Connect this product to only NLX BF (NLX nano motor), iMD (NBX motor) and Varios 170 (OPTION).
- When operating this product always consider the safety of the patient.
- Prior to use, always check for vibration, noise, and overheating. If any abnormalities are detected, stop using immediately and contact your Authorized NSK Dealer.
- Should the Product function abnormally, cease operation immediately and return the Product to your Authorized NSK Dealer for repair.
- Care should be taken not to place the motor cord near a gas burner. Never attempt to repair a burned motor cord. Always replace it with a new cord.
- Do not exceed the drive motor speed recommended by the manufacturer.
- In case of using a motor, ensure the Gear ratio on the display matches with the Gear ratio of the handpiece. Confirm that rotation speed is set in stipulated range.
- Do not use or leave the product in a high-temperature environment such as under strong direct sunlight, by a fire or near stove. It may cause a malfunction for inner circuit or sudden heat generation.
- Do not allow any impact on to the product. Do not drop the product.
- If chemical, solvent or antiseptic solution is adhered to this product, immediately wipe it away. Discoloration or deformation may occur if left as it is.
- Multi Pad is not suitable for autoclave sterilization. Do not autoclave Multi Pad nor sterilize it in any high temperature sterilization methods.
- If the unit smokes or smells of burning, turn off the power immediately and disconnect the power plug. Contact your Authorized NSK Dealer.

- The use of ACCESSORIES, transducers and cables, with the exception of transducers and cables sold by the manufacturer of this product as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.
- This product should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, this product should be observed to verify normal operation in the configuration in which it will be used.
- If any water drops remain on the handpiece or handpiece cord after autoclaving, wipe them off. Staining may result if left.
- Users are responsible for applying this product to a patient.
- Multi Pad needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information.
- Portable and mobile RF communications equipment can affect Electrical Medical equipment. Do not use RF equipment in close proximity to the product.
- LED lighting is built into the NBX motor and NLX nano motor. It may cause the poor visibility temporarily if you continue to look at the direct light of the LED. Please do not look at the LED light directly.

NOTICE

- Perform periodical maintenance checks.
- If the product has not been used for a long period, check for noise, vibration, and overheating before use.
- Responsibility for operating and maintaining medical devices belongs to the user.
- NLX BF/IMD are for cutting / polishing demanded by the general treatment of the tooth.
Varios 170 generates ultrasonic waves intended for use in dental applications such as scaling, root canal treatment, periodontal and cavity preparation.
- In case of repair, contact your Authorized NSK Dealer.

◆ Classifications of equipment

- Type of protection against electric shock:
 - Class II equipment
- Degree of protection against electric shock:
 - Type B applied part ⚠
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:
 - EQUIPMENT not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Mode of operation:
 - Continuous operation

Summary of user setting menu

• Motor

The following menus will be displayed on the LCD screen.

1. Language	▷	English
2. Gear ratio		Español
3. Lamp		Italiano
4. Endo mode		Français
5. Air		Deutsch
6. Endo alert		
7. Others		
8. Version		

1. Language			
2. Gear ratio	▷	Gear ratio 1	▷ 1 – 20 : 1 – 5 *
3. Lamp		Gear ratio 2	▷ 1 – 20 : 1 – 5 *
4. Endo mode			* Adjustable range
5. Air			
6. Endo alert			
7. Others			
8. Version			

1. Language				
2. Gear ratio				
3. Lamp	▷	On / Off	▷	On / Off
4. Endo mode		Intensity	▷	1.6 – 3.6 V * (NLX)
5. Air		Delay timer	▷	1.0 – 5.0 sec *
6. Endo alert				* Adjustable range
7. Others				
8. Version				

In case of iMD → 3.0 - 3.6 V * (iMD)

1. Language			
2. Gear ratio			
3. Lamp			
4. Endo mode	▷	Display setting	▷ Ncm
5. Air			Nmm
6. Endo alert			%
7. Others		Auto REV time	▷ 0.3 – 1.0 sec *
8. Version		Auto FWD time	▷ 1.0 – 3.0 sec
			* Adjustable range

"Endo mode" setting is not reflected in the case of the NLX nano motor.

1. Language			
2. Gear ratio			
3. Lamp			
4. Endo mode			
5. Air	▷	Display setting	▷ psi
6. Endo alert			bar
7. Others		Motor MIN pressure	▷ 0.03 – 0.10 MPa *
8. Version		Motor MAX pressure	▷ 0.20 – 0.40 MPa *
			* Adjustable range

1.Language
2.Gear ratio
3.Lamp
4.Endo mode
5.Air
6.Endo alert
7.Others
8.Version

▷ Endo mode alert	▷ On / Off
Alert setting 1	▷ 40 – 80 % *
Alert setting 2	▷ 80 – 100 % *
REV alert	▷ On / Off

* Adjustable range

*"Endo mode" setting is not reflected in the case of the NLX nano motor.

1.Language
2.Gear ratio
3.Lamp
4.Endo mode
5.Air
6.Endo alert
7.Others
8.Version

▷ LCD contrast	▷ -2, -1, 0, +1, +2
Backlight timer	▷ 1 – 30 min or ∞ *
Acceleration time (NLX)	▷ 0.5 – 3.0 sec *
Default setting	▷ No / Yes

In case of iMD

→ Cable length (iMD)	▷ 1.2 - 2.6m *
----------------------	----------------

* Adjustable range

1.Language
2.Gear ratio
3.Lamp
4.Endo mode
5.Air
6.Endo alert
7.Others
8.Version

▷ LCD XXXX	In case of iMD
NLX XXXX	→ iMD XXXX
VA170 XXXX	

• Scaler

In the case of connecting to the Varios 170(OPTION), this menu is displayed at the setup.

1.Version

▷ LCD XXXX	
VA170 XXXX	In case of iMD
NLX XXXX	→ iMD XXXX


1 FEATURES

- Electric Micromotor Integrated System (NLX BF / iMD) and Multi Function Ultrasonic Scaler (Varios 170) can be controlled with the Multi Pad system.
- User settings are achieved via the LCD graphical user interface.
- Visual identification of graphics on the Multi Pad display is possible even during operation of the handpiece.
- The Multi Pad allows up to 8 custom programs for your exacting needs (motor only).
- Multi Pad supports 5 languages (English, Español, Italiano, Français and Deutsch).

2 SPECIFICATIONS

Model	Multi Pad
Rated input	DC 12 V 100 mA 1.2 W
Dimensions	W95 x D138 x H31 mm
Weight	286 g (including cable and bracket)
Use Environment	Temperature : 10 – 40 °C Humidity : 30 – 75 %
Transportation and Store Environment	Temperature : – 10 – 60 °C Humidity : 10 – 85 % Atmospheric pressure : 500 – 1060 hPa

3 NAME OF EACH PART

4 START UP

When turning the Power On, the name of connecting instruments will be displayed for about 3 seconds along with a beep tone. Then, the LCD display will switch to operation screen automatically.

- NOTICE**
- When turning the Power On with the foot control being pushed, the display of the Multi Pad will show the error message with “Release Foot Pedal” along with beep tone. The error messages will disappear when the foot control is released.

5 DESCRIPTION FOR DISPLAY

5-1 MOTOR

Fig.1 Normal display when using micromotor

While the display shows the Fig. 1, you can control a motor by pressing the foot control pedal.

◆ Rotation display (Fig. 1-①)

When the motor is stopping, the display shows the Maximum Rotation Speed. When the Motor is rotating, the display shows the actual rotation speed. When the Motor stops, the display shows "0" on the display then turn to the Maximum Rotation Speed originally selected.

Forward Mode : Moves clockwise.

Reverse Mode : Moves counter clockwise.


◆ Rotation Speed display (Fig. 1-②)

The display shows Rotation Speed you selected. While the Motor is rotating, the actual rotation speed will be shown on the display.

◆ SET VALUE display (Fig. 1-③)

The pre-programmed Maximum Rotation Speed is indicated on the display.

NOTICE

-  icon will go off while the Motor is rotating. When increasing or reducing the Maximum Rotation Speed while the motor is rotating,  icon appears for 1 second.

◆ Program Number display (Fig. 1-④)

 The Program Number you selected will be displayed.
4 For detailed Program setting, refer to "7. USEFUL FUNCTIONS".

◆ Gear ratio display (Fig. 1-⑤)

 The preset Gear ratio is displayed. A total of 10 different Gear ratios can be selected. 8 Gear ratios are default and cannot be changed. 2 Gear ratios are for customized setting by users.
1 : 4 For detailed Gear ratio setting, refer to "6. SETTINGS".

5-2 SCALER

Fig.2 Normal display when using scaler

While the display shows the Fig. 2, you can control an ultrasonic scaler by pressing the foot control pedal.

◆ Ultrasonic Vibration display (Fig. 2- ①)

 When the scaler is vibrating, vibration icon will be displayed.

◆ Vibration Mode display (Fig. 2- ②)

Endo Selected Ultrasonic Vibration Mode is displayed. (Perio, Endo or General)

◆ Power Level display (Fig. 2- ③)

POWER 7  Selected Ultrasonic Vibration Power Level is displayed. (POWER 0 – 10)

5-3 SELECT THE INSTRUMENT

If 2 instruments (motor and scaler) are connected in a system, you can select the instrument with the following steps.

- 1) Ensure the Motor and Scaler have completely stopped.
- 2) Press **[ESC]** key at the normal display. The display will go to instrument selection mode when using Motor A or Scaler. The screen shown in Fig.3 is displayed.
- 3) Use **▲/▼** key to select the instrument you would like to use.
- 4) Press **[ENTER/SAVE]** key, then the instrument is selected.

NLX nano motor and Varios 170 (Option) are connected


Fig. 3 Display when selecting the instrument

6 SETTINGS

6-1 MOTOR

Functions you can set while the Motor is rotating: Maximum Rotation Speed, Rotation Direction.

Functions you can set while the Motor is stopping: Maximum Rotation Speed, Rotation Direction, Gear ratio, PRG.

◆ Speed setting

Press ▲/▼ key to set your desired speed at the normal display.

Speed Range: 1,000 – 40,000 min⁻¹.

NOTICE • Speed display changes faster as you press and hold the ▲/▼ key.

◆ Gear ratio setting

A maximum of 10 different Gear ratios can be set. 8 are fixed default settings, 2 for free settings.

- 1) Press **[SELECT]** key until 'Gear ratio' on the display blinks at the normal display.
- 2) Press ▲/▼ key to choose an appropriate Gear ratio.
- 3) Press and hold **[ENTER/SAVE]** key for 3 seconds or longer to save the set Gear ratio.

NOTICE • The turn selected by ▲/▼ key is shown below.

* Gear ratio 1 and 2 are only displayed when you set them in user setting menu.
For user setting menu, refer to "9. OTHER SETTINGS".

◆ Rotation Direction setting

Press **[FWD/REV]** key until "F" appears at the normal display for Motor Forward Rotation.

Press **[FWD/REV]** key until "R" appears at the normal display for Motor Reverse Rotation.

6-2 SCALER

◆ Treatment Mode setting

Press **[MODE]** key to select the different treatments at the normal display.

"Perio", "Endo" and "General" modes will be cycled through on every key press.

◆ Power Level setting

Press ▲/▼ key to select the Power Level at the normal display.

▲ key for increase, press ▼ key for decrease. (Power 0 to 10)

7 USEFUL FUNCTIONS

7-1 MOTOR

A maximum of 8 customized settings are available with the following selections.

- Maximum Rotation Speed
- Rotation Direction
- Gear ratio

◆ Program selection

- 1) Press **[SELECT]** key until program number displayed blinks at the normal display.
- 2) Press **▲/▼** key to select the program number that you want to use.
- 3) Press **[ENTER/SAVE]** key for 3 seconds or longer. The selection is completed.

NOTICE • When the program you selected is modified, the color of “PRG” will inverse.

8 ERROR CODES

8-1 NLX nano MOTOR

The system incorporates an automatic diagnosis feature that can help to diagnose the cause of the problem in the event of failure mode. When the system has failed, the display will show an “E-***” code and message. Details of each error code are listed below.

Error code	Error Message	Cause	Solution
E-00	Over current (Hard)	Excessive torque is loaded.	Remove the load from motor.
E-01	Over current (Soft1)	Abnormal electric current applied to motor and circuit.	Release the foot control.
E-02	Over current (Soft2)	Excessive current is applied for a given time.	Release the foot control.
E-03	Fault error	Overcurrent in the motor driver.	Contact your Authorized NSK Dealer.
E-04	Over heat	Motor is overused.	Leave the motor is until it cools down.
E-05	Input over voltage	Over voltage.	Contact your Authorized NSK Dealer.
E-06	Lamp over voltage	Lamp circuit over voltage.	Contact your Authorized NSK Dealer.
E-07	Residual voltage	Output circuit error.	Contact your Authorized NSK Dealer.
E-08	Over load	The load continued to exceed the limiter during the fixed time.	Remove the load from motor/handpiece, and release foot control.
E-09	Motor start failure	Handpiece/motor cord disconnection.	Check the handpiece cord is connected with the handpiece properly.
		Electric circuit may fail.	Contact your Authorized NSK Dealer.
E-10	Lamp under voltage	Lamp circuit under voltage.	Contact your Authorized NSK Dealer.
E-13	Esc error	Over motor control range.	Remove the load from motor/handpiece, and release the foot control.
E-14	EEPROM error	Memory read/write data error.	Contact your Authorized NSK Dealer.

8-2 NBX MOTOR

The system incorporates an automatic diagnosis feature that can help to diagnose the cause of the problem in the event of failure mode. When the system has failed, the display will show an “E-***” code and message. Details of each error code are listed below.

Error code	Error Message	Cause	Solution
E-00	Over Motor current	Excessive current is applied for a given time at high speed mode.	Remove the lead from motor/handpiece, and release foot control.
E-01	Over PAM voltage	Error in detecting the motor voltage more than DC28V.	Release the foot control.
E-02	PAM diff. failure	High speed Mode, Over difference voltage of PAM control.	Release the foot control.
E-03	Over EMF voltage	Low speed Mode, Over difference voltage of EMF control.	Release the foot control.
E-04	Torque diff. failure	Over difference voltage of Torque control.	Contact your Authorized NSK Dealer.
E-05	Input voltage error	Input voltage has exceeded the upper and lower limit values.	Check input voltage.
E-06	Over LED voltage	Detected over voltage of LED lamp.	Contact your Authorized NSK Dealer.
E-07	LED diff. failure	Over difference voltage of LED lamp control.	Contact your Authorized NSK Dealer.
E-08	Over load	Motor is overused.	Leave the motor until it cools down.
E-09	Motor start failure	Handpiece/motor cord disconnection.	Check the handpiece cord is connected with the handpiece properly.
E-11	Air FS failure	Air pressure is greater than 0.45[MPa].	Check user air pressure.
E-12	Comm. error	Disconnection or defective contact of COM port.	Check the connection or wiring of COM port.
E-13	Over current(endo)	Excessive current is applied for a given time at low speed mode.	Remove the lead from motor/handpiece, and release the foot control.
E-14	EEPROM error	Memory read/write date error.	Contact your Authorized NSK Dealer.

8-3 SCALER


Error code	Error Message	Cause	Solution
E-09	poorTool	Handpiece cord disconnection.	Check the handpiece cord is connected to the handpiece properly.
E-14	Fatal	Program error.	Contact your Authorized NSK Dealer.

9 OTHER SETTINGS

For each setting, follow the process below.

- 1) Ensure the Motor and the Scaler have completely stopped.
- 2) For normal display, press and hold **[SET UP]** key for 3 seconds or longer, then "User setting menu" will be displayed.
- 3) Use **▲/▼** key to select the item you wish to change.
- 4) Press **[ENTER/SAVE]** key, then select the item (Perform this process to change the settings after a setup indicated in 9-1 and 9-2").
- 5) Press and hold **[ENTER/SAVE]** key for 3 seconds or longer to save the program.

NOTICE • At the setting mode, when you move back to the previous screen, use the **[ESC]** key. If you press **[ESC]** key without setting, the menu below will appear. Follow the instruction on the display.

9-1 MOTOR SETTINGS

9-1-1 [1. Language] : Setting the display language.

1.Language	▷	English
2.Gear ratio		Español
3.Lamp		Italiano
4.Endo mode		Français
5.Air		Deutsch
6.Endo alert		
7.Others		
8.Version		

9-1-2 [2. Gear ratio] : 2 of 10 Gear ratio settings can be set for customized settings.

To move the cursor, press the **[SELECT]** key.

The Gear ratio can be set [20 to 1:1] or [1: 1 to 5].

The right side or the left side should be 1.


1.Language			
2.Gear ratio	▷	Gear ratio 1	▷ [1 - 20 : 1] or [1 : 1 - 5]*
3.Lamp		Gear ratio 2	▷ [1 - 20 : 1] or [1 : 1 - 5]*
4.Endo mode			
5.Air			
6.Endo alert			
7.Others			
8.Version			

* Adjustable range

9-1-3 [3. Lamp] : For setting the function of the LED.

At setup mode, select "Lamp" on the screen. Then you will see the next selection shown below diagram.


- On/Off : Setting the Motor Lamp on/off.
- Intensity : Setting the intensity of the Motor Lamp. 1.6-3.6V(NLX) / 3.0-3.6V(iMD). Defaults are 3.5V(NLX) / 3.3V(iMD).
- Delay timer : Setting the duration of lighting after the Motor is used (1.0 – 5.0 seconds). Default is 3.0 sec.

NOTICE • The setting items of "Intensity" and "Delay timer" will be displayed only when lamp setting is "On".

9-1-4 [4. Endo mode] : Endo (endodontics) mode setting.

This setting is not applicable to the NLX nano motor.


9-1-5 [5. Air] : Calibration for air driven foot control.

Air pressure may vary on the dental chair unit. The foot control of the Multi Pad system has a feature for synchronizing the motor speed with the air pressure.


Adjust the air pressure within a range shown below by adjusting foot control. When a message of "under" is displayed, air pressure is too low. When a message of "over" is displayed, air pressure is too high.

- Display setting : Changing air pressure unit display (psi, bar, MPa). Select the unit that you want to use.
- Motor MIN pressure : Setting the air pressure for the motor start-up.
(4.3 – 14.5 psi or 0.30 – 1.00 bar, 0.030 – 0.100MPa).
- Motor MAX pressure : Setting the air pressure for reaching the max Motor Rotation.
(29.0 – 58.0 psi or 2.00 – 4.00 bar, 0.200 – 0.400MPa).

9-1-6 [6. Endo alert] : Setting the alarm in Endo (Endodontics) mode.


These settings are not reflected in the case of the NLX nano motor.

9-1-7 [7. Others] : Other settings.


At setup mode, select “Others” on the screen. Then you will see the next selection shown below.

- LCD contrast : Setting the contrast of the LCD display.
Select the LCD contrast that you want to use, “- 2”, “- 1”, “0”, “+ 1”, “+ 2”. Default is “0”.
- Backlight timer : Setting the backlight timer of the LCD display.
Select the Backlight timer that you want to use, 1 minute to 30 minutes or infinite. Default is 10 min.
- Acceleration time : Setting the time for the motor rotation speed to reach the maximum speed value.
Select the acceleration time that you want to use, 0.5 seconds to 3.0 seconds. Default is 0.5 sec.
- Default setting : Back NLX BF / iMD to the factory default setting.
When you select “Default setting” by pressing **ENTER/SAVE** key, you will see the menu of Reconfirmation screen. Use **▲/▼** key to move to “Yes”, then **ENTER/SAVE** key for 3 seconds. Then the all settings are initialized.
- Cable length:Setting the cable length of NBX motor cord.
Select cable length from 1.2m to 2.6m. Default is 2.2m.

9-1-8 [8. Version] : Display product information.

- LCD : Software Version
 NLX/iMD : Software Version
 Varios 170(OPTION) : Software Version

9–2 SCALER SETTINGS

9–2–1 [1. Version] : Display product information.

LCD : Software Version

VA170 : Software Version

NLX/iMD : Software Version

10 DISPOSING PRODUCT

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

11 WARRANTY

Manufacture warrants this product to the Original purchaser against defects in material and workmanship under normal practices of installation, use and servicing. Article of consumption (such as O–Rings) is not covered by this warranty.

Symbols

Type B applied part.

Refer to Operation Manual.

Manufacturer.

Follow the waste of electric and electronic equipment (WEEE) Directive (2012/19/EU) to dispose of the product and accessories.

This conforms to CE European Directive of "Medical equipment directive 93/42/EEC."

Protected against vertically falling water drops.

Guidance and manufacturer's declaration – electromagnetic emissions

The Multi Pad is intended for use in the electromagnetic environment specified below. The customer or the user of the Multi Pad should assure that is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions EN 55011 CISPR 11	Group 1	The Multi Pad uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions EN 55011 CISPR 11	Class B	The Multi Pad is suitable for use in all establishments, including domestic establishments and those directly connected to the public low–voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions EN 61000-3-2 IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions EN 61000-3-3 IEC 61000-3-3	Complies	

EMC Measurement Conditions			
No.	Interface(s)	Max. cable length, shielding	Cable classifications
1	AC input (with plug)	Un-Specified	AC Power Supply Line
2	AC output (CN400)	4.0m, Shielded	AC Power Supply Line
3	DC power line (CN401 – CN105)	1.0m, Un-shielded	DC Power Supply Line
4	Motor line (CN104 – CN500)	1.0m, Un-shielded	AC Power Line
5	External I/F (CN300)	1.0m, Un-shielded	Input/Output Line
6	Control line for Multi Pad (CN302)	1.0m, Un-shielded	Input/Output Line
7	Lamp line (CN305 – CN501)	1.0m, Un-shielded	DC Power Line
8	Motor select line (CN306 – CN502)	1.0m, Un-shielded	DC Power Line
9	Motor line A (CN503)	2.2m, Un-shielded	AC Power Line
10	Lamp line A (CN5040)	2.2m, Un-shielded	DC Power Line
11	Motor line B (CN505)	2.2m, Un-shielded	AC Power Line
12	Lamp line B (CN506)	2.2m, Un-shielded	DC Power Line

Guidance and manufacturer's declaration – electromagnetic immunity

The Multi Pad is intended for use in the electromagnetic environment specified below. The customer or the user of the Multi Pad should assure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN 61000-4-2 IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst EN 61000-4-4 IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5 IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line to line ± 2 kV lines to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11 IEC 61000-4-11	< 5 % Ut (> 95 % dip in Ut) for 0.5 cycle 40 % Ut (60 % dip in Ut) for 5 cycles 70 % Ut (30 % dip in Ut) for 25 cycles < 5 % Ut (> 95 % dip in Ut) for 5 sec	< 5 % Ut (> 95 % dip in Ut) for 0.5 cycle 40 % Ut (60 % dip in Ut) for 5 cycles 70 % Ut (30 % dip in Ut) for 25 cycles < 5 % Ut (> 95 % dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Multi Pad requires continued operation during power mains interruptions, it is recommended that the Multi Pad be powered from an uninterruptible power supply or a battery.
Power frequency (50 / 60 Hz) magnetic field EN 61000-4-8 IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE : Ut is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The Multi Pad is intended for use in the electromagnetic environment specified below. The customer or the user of the Multi Pad should assure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			<p>Portable and mobile RF communications equipment should be used no closer to any part of the Multi Pad, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance. Where \sqrt{P} is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Conducted RF EN 61000-4-6 IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = 1.2\sqrt{P}$
Radiated RF EN 61000-4-3 IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz

NOTE 1 At 80MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobiles radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Multi Pad is used exceeds the applicable RF compliance level above, the Multi Pad should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Multi Pad.

b
Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Multi Pad

The Multi Pad is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Multi Pad can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of the transmitter W	Separation distance according to frequency of transmitter		
	150kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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NSK United Kingdom Ltd
NSK America Corp
NSK Oceania Pty Ltd