

POWER
SAFETY
ACCURACY
NO COMPROMISE.



Surgic Pro+

ORAL SURGERY
MICROMOTOR

POWERFUL
UP TO 80Ncm

**PLUS
3 FREE
Accessories**



POWER • SAFETY • ACCURACY

Surgic Pro+

MODEL	ORDER CODE	RRP PRICE
Surgic Pro Optic	Y1001933	\$9,549

INCLUDES

- Implant handpiece **TiMax X-SG20L**
- Micromotor **SGL70M**

**SAVE
\$2,750**



**FREE
Travel Case**

ORDER CODE Y1001952
RRP-\$831



**FREE
Sterilisation
Cassette**

MODEL **SG-CASE**
ORDER CODE Z313101
RRP-\$691



**FREE
Irrigation Tube**

BOX OF 50
ORDER CODE Y900113
RRP-\$1,228

Surgical Micromotor & Cable



**FREE
Sterilisation
Cassette**



Optic	ORDER CODE	RRP PRICE
SGL70M (with 2m Cord)	E1023001	\$4,885
SGL70M 3.5m (with 3.5m Cord)	E1023002	\$5,316
SG-CASE	Z313101	FREE RRP-\$691

Non-Optic	ORDER CODE	RRP PRICE
SG70M (with 2m Cord)	E1025001	\$3,883
SG70M 3.5m (with 3.5m Cord)	E1025002	\$4,313
SG-CASE	Z313101	FREE RRP-\$691

iCart Duo

The NSK iCart Duo allows you to install both the NSK SurgicPro and Variosurg Systems as well as additional accessories on to the one easily transportable cart.

MODEL	ORDER CODE	PROMO PRICE
iCart Duo	S9090	\$1,809 RRP-\$2,010

Specifications

- Dimensions : H101.65 cm
- Weight : 16.5 kg
- Adjustable Height
- Lockable Wheels

**10%
OFF**



"As my business grows NSK has always supported my practice.



My Surgic Pro is consistently reliable, created with precision, and a pleasure to work with.

It is the workhorse of my practice and I couldn't live without it!"

Dr Cecilia So BDS Hons (USyd), MFDS (Lond), Grad Dip Clinical Dentistry (Oral Implants)(USyd), FICD, Clinical Associate Lecturer (Oral Implants) (USyd)

Osseo 100 + Five FREE MultiPegs

OSSEOINTEGRATION MONITORING DEVICE

MODEL	ORDER CODE	RRP PRICE
Osseo 100	Y1004175	\$3,575




**SAVE
\$355**

NSK's Osseo 100 measures implant stability and osseointegration to enhance decisions about when to load the implant. This is especially important when using protocols with shorter treatment time and treating higher risk patients.

The uncomplicated operation that measures ISQ, allows the implant loading period to be planned in advance.

3-STEP PROCEDURE



The MultiPeg™ is attached to the implant. It screws effortlessly into the implant's internal threads (approximately 6-8 Ncm of torque).



Just aim for the magnet on top of the MultiPeg™. Non-invasive, objective, accurate and repeatable. The peg is excited by magnetic pulses and vibrate due to the stiffness in the contact area between the bone and the implant surface.



An ISQ value is generated and shown on the display. This reflects the level of stability on the universal ISQ scale – from 1 to 99. The higher the ISQ value, the more stable the implant.

MulTipeg™ REUSABLE

- MulTipeg™
- For all major implant systems *
- Tissue friendly, durable titanium
- Autoclavable appx. 20 times
- Optimal platform fit
- ISQ Standard Calibrated

*There are different MulTipegs™ available made to fit different implant system and types. Please refer to the MulTipeg list from the supplier, or scan the QR code to view online.

MODEL	ORDER CODE	PROMO PRICE
MulTipeg™	Various	FREE RRP \$71 EACH



Scan the QR code for multipeg list



Implant Handpieces



MODEL **Ti-Max X-SG20L**
ORDER CODE C1003

- Optic
- Reducing 20:1
- Titanium
- CA Burs

PROMO PRICE
\$2,381

RRP \$2,976



MODEL **Ti-Max X-DSG20L**
ORDER CODE C1068

- Optic
- Reducing 20:1
- Titanium
- Dismantable
- CA Burs

PROMO PRICE
\$2,567

RRP \$3,208



MODEL **S-Max SG20**
ORDER CODE C1010

- Non Optic
- Reducing 20:1
- Stainless Steel
- CA Burs

PROMO PRICE
\$1,451

RRP \$1,813



Speed Increasing 1:3



MODEL **Ti-Max X-SG93L**
ORDER CODE C1004

- Optic
- Increasing 1:3
- Triple spray
- FG Burs

PROMO PRICE
\$2,287

RRP \$2,859



MODEL **Ti-Max Z-SG45L**
ORDER CODE C1107

- Optic
- Increasing 1:3
- Triple Spray
- 25mm FG Burs

PROMO PRICE
\$1,944

RRP \$2,430

Direct Drive 1:1



MODEL **Ti-Max X-SG25L**
ORDER CODE C1011

- Optic
- Direct Drive 1:1
- CA Burs

PROMO PRICE
\$2,107

RRP \$2,633



MODEL **Ti-Max X-SG65L**
ORDER CODE H1009

- Optic
- Direct Drive 1:1
- CA Burs

PROMO PRICE
\$2,075

RRP \$2,593

Microsurgical Handpieces



MODEL **SGS-ES**
ORDER CODE H264

- Non Optic
- Direct Drive 1:1
- 65mm HP Burs

PROMO PRICE
\$1,331

RRP \$1,663



MODEL **SGA-E2S**
ORDER CODE H265

- Non Optic
- Increasing 1:2
- 65mm HP Burs

PROMO PRICE
\$2,164

RRP \$2,705



MODEL **SGA-ES**
ORDER CODE H263

- Non Optic
- Direct Drive 1:1
- 20° angulation
- 65mm HP Burs

PROMO PRICE
\$1,845

RRP \$2,306



MODEL **SGS-E2S**
ORDER CODE H266

- Non Optic
- Increasing 1:2
- 20° angulation
- 65mm HP Burs

PROMO PRICE
\$1,560

RRP \$1,950

ISD900

Faster and Safer Implant Treatments

NSK's iSD900 cordless screwdriver helps to safely place and remove cover screws, healing caps and abutments during dental implant procedures, making treatment up to 50% faster.

MODEL	ORDER CODE	PROMO PRICE
ISD900	Y1001358	\$3,154
		RRP \$3,943

- Universal to major implant systems
- Wide torque range 10-40 Ncm
- Adjustable in 1 or 5 Ncm increments
- Speed range 15, 20, 25 min⁻¹



iSD900


DualDefense™ Proprietary Infection Control by NSK

Clean Head System

Anti-suck back system prevent cross contamination of the air lines

Anti-Retraction Valve

Anti-Retraction systems prevents cross contamination of the water lines in High-Speed Air Turbines and Electric Micromotors.



Ti-Max Z LED Premium Titanium 2 YEAR WARRANTY

20%
OFF



MODEL **Ti-Max Z900L**
ORDER CODE P1111
• 26W Power
• Standard Head
PROMO PRICE
\$1,735
RRP \$2,168

MODEL **Ti-Max Z800L**
ORDER CODE P1110
• 23W Power
• Miniature Head
PROMO PRICE
\$1,735
RRP \$2,168

S-Max M Advanced Stainless Steel 1 YEAR WARRANTY



MODEL **S-Max M900L**
ORDER CODE P1254
• 26W Power
• Standard Head
PROMO PRICE
\$1,210
RRP \$1,512

MODEL **S-Max M800L**
ORDER CODE P1253
• 23W Power
• Miniature Head
PROMO PRICE
\$1,210
RRP \$1,512



S-Max Pico Paediatric Highspeed ULTRA MINI HEAD 1 YEAR WARRANTY

MODEL **S-Max Pico**
ORDER CODE P1140
• 8W Power
• Ultra Short Burs available
PROMO PRICE
\$1,210
RRP \$1,512

Non-Optic options available.
Other connections offered for KaVo, Sirona, W&H, Bien-Air.



Ti-Max X450L Surgical Highspeed REAR EXHAUSTED AIR 2 YEAR WARRANTY

MODEL **Ti-Max X450L**
ORDER CODE P1078
• 21W Power
• 45° head
• Triple Water Jets
PROMO PRICE
\$1,468
RRP \$1,835

Buy ANY 6 Ti Max Handpieces Get FREE iCare*

*Buy any 6 Ti Max handpieces at RRP and receive a FREE iCare and adaptors of your choice!

- Pre-set solution volume based on handpiece type
- Thorough cleaning and lubrication with simple operation
- Chuck cleaning and lubricating function
- Purging mode removes excess oil
- Integrated 1.2 litre oil tank reduces burden of refill



iCare



iCare C2 Type

ORDER CODE Y1002787

Dual high speed + Dual low speed Connectors requires 2 high speed Adaptors*



iCare C3 Type

ORDER CODE Y1002788

Single high speed + Triple low speed Connectors requires 1 high speed Adaptor*



*Adaptors available to fit NSK, KaVo, Sirona, Bien-Air, Morita high speed handpieces. Limit x1 free adaptor free per each iCare C3 and x2 free adaptors per each iCare C2.

Varios Combi Pro 2 Maintenance Systems in 1

Buy Varios Combi Pro Basic set and get 2 more Prophy handpieces, another scaler handpiece and 4 bottles of powder FREE!

Basic Set

MODEL Varios Combi Pro Basic Set
ORDER CODE Y1002843

CONTENTS

- Prophy-Handpiece with 60° Nozzle
- Prophy-Hose
- Prophy-Chamber
- FLASH pearl (5pcs. x 15g sachets)
- Maintenance kit
- VA2Lux-Handpiece
- VA2Lux-Hose
- 3 Scaling tips
- Sterilization case
- Irrigation bottle (400mL)
- Foot control



MODEL	ORDER CODE	RRP PRICE
Varios Combi Pro Basic Kit	Y1004194	\$7,274
Prophy HP x 2	T1074001	FREE RRP-\$398-EACH
Varios2 Lux Handpiece x 1	E351050	FREE RRP-\$1,905
Flash Pearl Powder Box of 4 bottles	Y900698	FREE RRP-\$205-EACH

**\$7,274
with FREE
Handpieces
& Powder**



NSK Oceania Ltd www.nskoceania.co.nz

Product Information | Technical Support
Freephone: 0800 44 33 21 (New Zealand) Email: info@nskoceania.com.au