

THE NEW STYLE ELECTRIC THE COMPACT MAKES A BIG CHANGE

NLZ E ELECTRIC MICROMOTOR SYSTEM



NEW PRODUCT

A totally new installation style only from NSK



Specification comparison

| | | Our existing model | NLZ • NLZ E |
|--|---------------------|--------------------------------|--|
| Motor specs | Number of rotations | 2,000-40,000 min ⁻¹ | 100-40,000 min⁻¹ |
| | Max torque | 3.4 Ncm | 4.2 Ncm 120% High Torque |
| Corresponding model to endodontic mode | | - | ○ (NLZ E) |

NLZ-E

Electric Micromotor

PORTABLE SYSTEM WITH ENDO FUNCTION

+ FREE NLZ Endo Handpiece and ISO-E Lux Probe

- Easily convert air driven handpieces to electric with the all new NLZ-E.
- Correctly transmit reciprocating actions for a smooth endodontic treatment.
- New NLZ ENDO 6:1 reciprocating handpiece.
- Endo torque calibration function – This function measures the 'no-load' current values and corrects the torque accordingly, resulting in greatly enhance safety in endodontic treatment.
- Contra-Check Function helps detect abnormalities caused by deterioration of the contra-angle handpiece and insufficient maintenance before use.

| MODEL | ORDER CODE | RRP PRICE |
|--------------------|------------|--------------------------|
| NLZ-E Portable | Y1500361 | \$6,502 |
| NLZ Endo Handpiece | C1130 | FREE RRP-\$977 |
| ISO E-Lux Probe | Z1265 | FREE RRP-\$814 |



NLZ E





NLX nano

Electric Micromotors

PORTABLE SYSTEM WITH NANO CONTROL

This system can help any user of air driven handpieces to convert to an electric system.




| MODEL | ORDER CODE | PROMO PRICE |
|--------------|------------|-------------------------------|
| NLX Portable | Y1001362 | \$4,859 RRP \$6,074 |



INTEGRATED SYSTEM WITH MULTIPAD CONTROL

This system can be easily integrated into your existing chair to control up to 2 electric micromotors and 1 Varios scaling system.

| MODEL | ORDER CODE | PROMO PRICE |
|---------------------------------|----------------------|-------------------------------|
| NLX nano standard with MultiPad | NLXNANOIntegratedpad | \$4,145 RRP \$5,181 |
| NLX plus with MultiPad | NLXPLUSIntegratedpad | \$4,363 RRP \$6,454 |



Varios 170 iPiezo engine



Varios 170

LED Illumination

INTEGRATED MULTI FUNCTION ULTRASONIC SCALER

The smallest module for easy integration into new and existing dental equipment. Contents: Built in module, Handpiece, Handpiece cord, 3 scaling tips (G4, G6, G8), 3 x Tip replacement wrench with torque limiter, Sterilisation case. Installation not included.

| MODEL | ORDER CODE | PROMO PRICE |
|-----------------------------------|------------|-------------------------------|
| Varios 170 Optic Complete Set | Y1001351 | \$1,977 RRP \$2,472 |
| Varios 170 Non-Optic Complete Set | Y1001352 | \$1,404 RRP \$1,755 |

20% OFF

S-Max M Handpieces

Advanced Stainless Steel



MODEL **S-Max M900L**
ORDER CODE P1254
• 26W Power
• Standard Head
PROMO PRICE \$1,210
RRP \$1,512



MODEL **S-Max M800L**
ORDER CODE P1253
• 23W Power
• Miniature Head
PROMO PRICE \$1,210
RRP \$1,512



MODEL **S-Max M65**
ORDER CODE H1008
• Non-Optic
• 1:1 Direct Drive
• Single spray
• HP/CA burs
PROMO PRICE \$716
RRP \$895



Mini Air Micro Motors

MODEL **S-Max M205**
ORDER CODE M1007001
• Non-Optic
• 1:1 Direct Drive
• Internal Spray
• Non retraction valve
PROMO PRICE \$1,063
RRP \$1,328

Non-optic options available. Other connections offered for KaVo, Sirona, W&H, Bien Air.



MODEL **S-Max M95L**
ORDER CODE C1023
• Optic
• Increasing 1:5
• Quattro spray
• FG Burs
PROMO PRICE \$1,543
RRP \$1,929



MODEL **S-Max M25L**
ORDER CODE C1024
• Optic
• 1:1 Direct Drive
• Single spray
• CA burs
PROMO PRICE \$1,184
RRP \$1,480



MODEL **S-Max M15L**
ORDER CODE C1025
• Optic
• 4:1 Reduction
• Single Spray
• CA Burs
PROMO PRICE \$1,441
RRP \$1,802

"The Osseo 100 unit has been extremely easy and affordable to use. Being able to compare RFA readings gives me confidence that we are ready to progress to the restorative phase."



Dr David Attia BOH (DentSci) (Griffith), GradDipDent (Griffith), MSc (Oral Implantology) (Goethe), PGDipClin (Orth)(CoL), FICOI

Osseo 100 + FREE Sterile covers

Osseo 100 measures implant stability and osseointegration using RFA, to enhance decisions about when to load the implant.

MODEL **Osseo 100**
ORDER CODE Y1004175
RRP PRICE \$3,575

MODEL **Osseo 100 Sterile Covers**
ORDER CODE 55105
PROMO PRICE FREE
RRP \$417

FREE Sterile Covers THREE PACKS OF 20



Simple and precise root canal treatment

Endodontic Motors

ENDO-MATE TC2

MODEL **ENDO-MATE TC2 (230V)**
MP-F16R head
ORDER CODE Y1001027

- Large LCD display with 5 programs for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1 to 4.0 NCM depending on head:gear ratio

PROMO PRICE \$1,981
RRP \$2,476



ENDO-MATE AT

MODEL **Endo-Mate AT (230V)**
ORDER CODE Y1002190

- Fits all major brands of Files
- Torque control, with integrated auto-reverse function
- Compact stylish control/monitor pad for simple and user friendly operation

PROMO PRICE \$1,583
RRP \$1,979



Apex Locator

IPEX II

MODEL **iPEX II Complete set**
ORDER CODE Y1002208

- Accurate measurement of any root canal
- SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals

PROMO PRICE \$1,484
RRP \$1,855



20% OFF

"The NSK Endomate TC2 has been my motor of choice for the past 10 years. It is a well balanced, cordless motor that delivers consistent torque and speed for the endodontic treatment."



Dr Morgan Wu
Specialist Endodontist
Bds (Otago), Dcd (Endo)(Melb)

Varios Combi Pro Basic Kit

+ FREE Perio Kit

2 Maintenance Systems in 1

Varios Combi Pro combines ultrasonic scaling and powder therapy for a more effective and comfortable oral hygiene and biofilm control

Varios Combi Pro

Basic Kit

MODEL Varios Combi Pro Basic Kit
ORDER CODE Y1002843

- Advanced optic iPiezo scaling for General, Perio & Endo treatment.
- Powder Therapy to clean Supragingival areas efficiently and quickly.
- User friendly handpiece and nozzle for easy access.

Perio Kit

MODEL Varios Combi Pro Perio Kit
ORDER CODE Y1003042

- Remove biofilm from Subgingival areas with the gentle glycine based Perio-Mate Powder.
- Sterile single-use plastic nozzle tips.
- Safely remove biofilm inside periodontal pockets and around dental implants.



FREE PERIO KIT WORTH \$2,748

| MODEL | ORDER CODE | RRP PRICE |
|----------------------------|------------|----------------------------|
| Varios Combi Pro Basic Kit | Y1004194 | \$7,604 |
| Varios Combi Pro Perio Kit | Y1003042 | FREE RRP \$2,748 |



Scan the QR code for more information

Varios2 Scaler Handpiece

| MODEL | ORDER CODE | PROMO PRICE |
|---------------------------------------|------------|-------------------------------|
| Varios2 Lux HP Optic Scaler Handpiece | E351050 | \$1,619 RRP \$1,905 |
| Varios2 HP Non-Optic Scaler Handpiece | E350050 | \$987 RRP \$1,161 |



15% OFF VARIOS HANDPIECES

Wide choice of Ultrasonic Scaler Tips

Buy 4 tips + FREE torque wrench

| MODEL | RRP PRICE |
|------------------------|-------------------|
| Ultrasonic Scaler Tips | From \$124 |

| MODEL | ORDER CODE | PROMO PRICE |
|----------------------|------------|--------------------------|
| Varios Torque Wrench | Z221076 | FREE RRP \$146 |



Scan the QR code for full tip guide



BUY 4 TIPS GET FREE TORQUE WRENCH

MOBILE DENTISTRY

VIVA ace

+ FREE ISO E-Luxe Probe

MOBILE DENTISTRY SYSTEM

For more comfortable and more advanced home visits for dental treatment. VIVA ace combines the functionality of a dental treatment centre in a compact body with multipurpose micromotor, ultrasonic scaler and more. Its lightweight build means it is easy to carry and straight forward to set up in different locations. The first of its kind with high-functionality, supports more comfortable dental treatment wherever the patient.

Price excludes handpieces.

| MODEL | ORDER CODE | PROMO PRICE |
|-----------------------|------------|---------------------------------|
| VIVA ace complete set | Y1004228 | \$19,800 RRP \$22,000 |



ISO E-Lux Probe

INTENSE DIRECT ILLUMINATION

Intense direct illumination for diagnosis of cracks in teeth plus occlusal and interproximal caries not visible to the naked eye.

| MODEL | ORDER CODE | PROMO PRICE |
|-----------------|------------|---------------------------|
| ISO E-Lux Probe | Z1265 | \$651 RRP \$814 |

FREE
with
VIVA ace
purchase

VIVA ace

MOBILE DENTISTRY SYSTEM

Ask about
our bundle
options



CORELESS MICROMOTOR

VOLVERE i7

A USER-FRIENDLY MICROMOTOR FOR ANY DENTAL LABORATORY

The VOLVERE i7 is another advance in NSK's quest to provide dental technicians access to its globally renowned dental laboratory engines. The stylish VOLVERE i7 delivers a compact, advanced design, all essential features and excellent value for money.

- Compact, smart, light, flexible, torque 4.1 NCM
- Fits into the smallest space
- Microcomputer Control
- Auto Cruise Function

20% OFF

| MODEL | ORDER CODE | PROMO PRICE |
|---|------------|-------------------------------|
| NSK Volvere i7 RM Laboratory Micromotor System 35,000 rpm complete with Standard Micromotor | Y1003493 | \$1,852 RRP \$2,315 |



VOLVERE i7



Standard Micromotor

NSK Oceania Ltd www.nskoceania.co.nz

Product Information | Technical Support
Freephone: 0800 44 33 21 (New Zealand) Email: info@nskoceania.com.au

All offers as listed. Valid until the 30th September 2021. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in NZD and exclude GST.