



Special Offers

APRIL - JUNE 2025

Smart and intuitive.

NEW

VarioSurg 4

The Essential Ultrasonic Bone Surgery System for Implant Treatments

VarioSurg 4 ULTRASONIC BONE SURGERY SYSTEM

NEW

BUY VARIOSURG 4 **GET FREE ADDITIONAL** OPTIC SCALER WITH 2M CABLE WORTH \$6,521

Perform safer bone surgery with NSK VarioSurg 4

With its stylish design and wireless foot control, VarioSurg 4 builds on the reliability and control of previous models to deliver a modern ultrasonic bone surgery system, with Piezo technology to minimise damage to surrounding tissue and bone.

Specifications – Control Unit

- Frequency: 28-32 kHz
- Power supply: AC100 - 240V 50/60 Hz Max.
- Pump output 76 mL/min
- Programs: SURG x5, ENDO x2, PERIO x2
- Dimensions: W 245 x D 235 x H 90 mm

Complete set

MODEL VARIOSURG 4 (230V)
ORDER CODE Y1500768

RRP \$14,735

Includes:

- Control unit
- LED handpiece with 2m cord
- FC-86 foot control
- Sterilisation cassette
- Handpiece stand
- Irrigation tube (3 pcs.) and other accessories
- Basic H-S Kit (H-SG1, SG3, SG5, SG6D, SG7D, SG11 and Tip holder)



Complete set without Foot Control

MODEL VARIOSURG 4 W/O FC (230V)
ORDER CODE Y1500769

RRP \$11,890

Includes:

- Control unit
- LED handpiece with 2m cord
- Sterilisation cassette
- Handpiece stand
- Irrigation tube (3 pcs.) and other accessories
- Basic H-S Kit (H-SG1, SG3, SG5, SG6D, SG7D, SG11 and Tip holder)

For more information scan the QR code.



Surgic Pro2 IMPLANT AND ORAL SURGERY

Specifications

- Power supply: AC 100-240V 50/60 Hz Max
- Pump output: 75 mL/min
- Programs: 8 Programs/Implant Systems
- Dimensions: W 245 x D 235 x H 90 mm
- Torque: 5-80 Ncm
- Motor speed: 200-40,000 min⁻¹
- Coolant solution flow volume button
- Light power (optic micromotor) over 32,000 LUX
- Foot control functions
- PRG (Program) button
- Forward/Reverse button
- Speed control pedal

Complete set with X-SG20L

MODEL SURGIC PRO2 OPT
ORDER CODE Y1500511

Includes:

- Control Unit
- Unit SGL80M Optic Micromotor
- X-SG20L Optic Handpiece (20:1 Reduction)
- Foot control wireless
- Irrigation tube (pack of 3) and other accessories

RRP \$10,313

PROMO PRICE \$8,250

SAVE \$2,063

For more information scan the QR code.



OSSEO 100+ **20% OFF**

Osseo 100+ OSSEOINTEGRATION MEASURING DEVICE

Removes doubt

The 'Osseo 100+' measures the implant stability quotient contactlessly and can be connected using Bluetooth®, permitting ISQ confirmation over a shared interface. You can share and manage the data of those measured ISQ scores on other terminals through Surgic Pro2. Multipegs sold separately.



MODEL OSSEO 100+
ORDER CODE Y1004176

RRP \$4,706

PROMO PRICE \$3,764

SAVE \$942

For more information scan the QR code.



Surgical Handpieces

Wide range of surgical handpieces, to suit a variety of procedures, from implants, to surgical extractions and more.

For more information scan the QR code.



20% OFF ALL SURGICAL HANDPIECES



Implant Handpiece
MODEL S-MAX SG20
ORDER CODE C1010

- Non-Optic
- Reducing 20:1
- Stainless Steel
- CA Burs

RRP \$1,959

PROMO PRICE \$1,567

SAVE \$392



Micro Surgery Handpiece
MODEL SGA-ES
ORDER CODE H263

- Non-Optic
- Direct Drive 1:1
- 20° angulation
- 65mm HP burs

RRP \$2,491

PROMO PRICE \$1,993

SAVE \$498



Micro Surgery Handpiece
MODEL SGS-ES
ORDER CODE H264

- Non-Optic
- Direct Drive 1:1
- 65mm HP Burs

RRP \$1,797

PROMO PRICE \$1,438

SAVE \$359



Rear Exhaust Handpiece
MODEL Ti-MAX X450L
ORDER CODE P1078

- Optic
- Rear-exhausted 45° head
- 21W Power

*Standard head for NSK coupling

RRP \$1,982

PROMO PRICE \$1,586

SAVE \$396



Rear Exhaust Handpiece
MODEL X450KL
ORDER CODE P1079

- Optic
- Titanium high speed HP
- 45° angled head

*Standard head for KaVo®.

RRP \$2,193

PROMO PRICE \$1,754

SAVE \$439



Contra-angle Handpiece
MODEL Ti-MAX Z45L
ORDER CODE C1064

- Optic
- 1:4.2 Speed Increasing
- Switch Spray function
- 45° angled head
- For FG burs



RRP \$2,675

PROMO PRICE \$2,140

SAVE \$535

High-Speed Handpieces

Ti-MAX Z PREMIUM

For more information scan the QR code.



20% OFF A SINGLE HANDPIECE – OR BUY 2 @ RRP AND GET 1 FREE (OF EQUAL OR LESSER VALUE)

Dynamic Power System

Equipped with a new system, the "DYNAMIC POWER SYSTEM" that provides a revolutionary cutting feel. The enlarged exhaust port, optimized nozzle structure, and enlarged rotor chamber achieve the highest rotational power in its category.

Specifications

- Titanium Body
- DURAGRIP
- Cellular Glass Optics
- Ceramic Bearings
- Clean Head System
- Push Button Chuck
- Quattro Spray

- Quick Stop
- Microfilter*
- DYNAMIC POWER SYSTEM*
- DURAPOWERC HUCK*
- Simple Cartridge Change

*Only Z990L / Z890L



Miniature Head Handpiece
MODEL Ti-MAX Z890L
ORDER CODE PA23720001

- Power: 41 W**
- Speed: 320,000-400,000 min⁻¹
- Head Size: ø10.6 x H 10.8 mm

**At 3-4.2 bar.

RRP \$2,342

PROMO PRICE \$1,874

SAVE \$468



Standard Head Handpiece
MODEL Z990L
ORDER CODE PA23730001

- Power: 44 W**
- Speed: 280,000-360,000 min⁻¹
- Head Size: ø12.1 x H 12.7 mm

**At 3-4.2 bar.

RRP \$2,342

PROMO PRICE \$1,874

SAVE \$468



Micro Head Handpiece
MODEL Ti-MAX Z MICRO
ORDER CODE PA2357

- Power: 20 W
- Speed: 380,000-450,000 min⁻¹
- Head Size: ø9.0 x H 10.8 mm

RRP \$1,633

PROMO PRICE \$1,306

SAVE \$327

Ti-Max Z890 Other Connections

- KaVo®, Dentsply Sirona® & W&H®

RRP \$2,452

PROMO PRICE \$1,962

SAVE \$490

Ti-Max Z990 Other Connections

- KaVo®, Dentsply Sirona® & W&H®

RRP \$2,452

PROMO PRICE \$1,962

SAVE \$490

Ti-Max Z Micro Other Connections

- KaVo®, Dentsply Sirona® & W&H®

RRP \$1,881

PROMO PRICE \$1,505

SAVE \$376

NLZ Electric Micromotor

WITH Ti-MAX Z95L

GET 25% OFF NLZ PORTABLE AND Ti-MAX Z95L

Switch to an electric motor with this easy plug and play system. Boost your crown preps and take control with the NLZ Motor and Ti-Max Z95L handpiece.

Brushless Electric Micromotor Portable System

MODEL **NSK NLZ S230**
ORDER CODE **Y1500360**

Specifications:

- Optic
- Speed: 2,000-40,000 min⁻¹
- Torque: 3.4 Ncm
- Internal water spray
- Rated input: AC 28V 50/60 Hz
- Dimensions: W127 x D149 x H54 (mm)
- Weight: 517 g

Includes:

- Control & Main unit
- AC adaptor
- NLZ Motor/cord/bracket

Does not include handpiece.

RRP \$5,618

PROMO PRICE \$4,213

SAVE \$1,405



Ti-Max Z Series

As specialists in ultra-high-speed rotary equipment, NSK is passionate about creating innovative products that deliver outstanding value. NSK manufactures almost all parts in-house to guarantee the highest precision, reliability and unmatched performance of all products.

The Ti-Max Z contra-angle and air turbine series resulted from NSK's constant quest to make the impossible possible.

MODEL **Ti-MAX Z95L**
ORDER CODE **C1034SW**

- Optic
- 1:5 Speed Increasing
- Switch Spray function
- For FG Burs

RRP \$2,507

PROMO PRICE \$1,880

SAVE \$627



For more information scan the QR code.



20% OFF SELECTED S-MAX RANGE

Stainless steel handpieces, all with internal water and choice of Optic or Non-Optic.

S-Max M Series

The S-Max M series is the next phase of evolution of our globally acclaimed handpieces. These new models are the embodiment of NSK's expertise in handpieces. The M series delivers the reliability and safety demanded by both practitioners and the ever-changing standards of the dental instrument industry.

MODEL **S-MAX M900L**
ORDER CODE **P1254**

- Optic
- 26W Power
- Stainless Steel
- Standard head

RRP \$1,633

PROMO PRICE \$1,306

SAVE \$327

MODEL **S-MAX M900**
ORDER CODE **P1256**

- Non-Optic
- 26W Power
- Stainless Steel
- Standard head

RRP \$1,545

PROMO PRICE \$1,236

SAVE \$309



MODEL **S-MAX M95L**
ORDER CODE **C1023**

- Optic
- 1:5 Increasing
- Quattro Water Spray
- Stainless Steel
- For FG burs (ø1.6)

RRP \$2,084

PROMO PRICE \$1,667

SAVE \$417

MODEL **S-MAX M95**
ORDER CODE **C1026**

- Non-Optic
- 1:5 Increasing
- Quattro Water Spray
- Stainless Steel
- For FG burs (ø1.6)

RRP \$1,661

PROMO PRICE \$1,329

SAVE \$332



MODEL **S-MAX M25L**
ORDER CODE **C1024**

- Optic
- 1:1 Direct Drive
- Single Spray
- Stainless Steel
- For CA burs

RRP \$1,599

PROMO PRICE \$1,279

SAVE \$390

MODEL **S-MAX M25**
ORDER CODE **C1027**

- Non-Optic
- 1:1 Direct Drive
- Single Spray
- Stainless Steel
- For CA burs

RRP \$1,122

PROMO PRICE \$898

SAVE \$224



E-Type Non-Optic Airmotor

MODEL **S-MAX M205**
ORDER CODE **M1007001**

- Non-Optic
- Internal Water Spray
- Max. 25,000 min⁻¹
- Midwest 4-hole port

RRP \$1,435

PROMO PRICE \$1,148

SAVE \$287



Connects to



For more information scan the QR code.



Perio-Mate

BIOFILM ERASER

15% OFF FOR PERIO-MATE AND PROPHY-MATE NEO COMPLETE SYSTEMS

Range of applications from supragingival to subgingival use

- Treatment of periodontal pockets and peri-implantitis
- Eradication of biofilm in periodontal pockets
- Removal of plaque and discoloration
- NSK's innovative nozzle tip
- Precise, flawless powder flow
- Slim nozzle design
- Smooth 360° rotation
- Powder control ring

MODEL PER-PTL-P for NSK FlexiQuik Coupling
ORDER CODE Y1002654

Includes:

- Handpiece
- Powder case
- Nozzle tip (40 pcs.)
- Maintenance kit
- Perio-Mate Powder 160g x 1

RRP \$2,005

PROMO PRICE \$1,704

SAVE \$301

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®.



For more information scan the QR code.



Prophy-Mate neo

AIR POLISHING SYSTEM

The best polishing tool for beautiful white teeth

- Removal of plaque and discoloration
- Preparation for pit and fissure sealant
- Simple to operate
- Easy maintenance
- Sturdy and reliable
- Quick connect to major brands of turbine couplings
- Comes with 60° and 80° nozzles

MODEL PMNG-PTL-P for NSK FlexiQuik Coupling
ORDER CODE Y135029

Includes:

- Handpiece with 60° nozzle head
- Powder case
- Replacement nozzle head 80°
- 5 pcs. x 15 g sachets of FLASH pearl cleaning powder
- Maintenance kit

RRP \$2,005

PROMO PRICE \$1,704

SAVE \$301

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®.



For more information scan the QR code.



Varios Combi Pro

VARIOS COMBI PRO BASIC SET & VARIOS COMBI PRO PERIO SET

15% OFF VARIOS COMBI PRO BASIC SET & PERIO SET

MODEL VARIOS COMBI PRO BASIC SET
ORDER CODE Y1004194

Includes:

- PROPHY-Handpiece with 60° Nozzle
- PROPHY-Hose
- PROPHY-CHAMBER
- FLASHpearl (5 pieces x 15g sachets)
- Maintenance kit
- VA2LUX-Handpiece
- VA2LUX-Hose
- 3 Scaling tips
- Sterilisation case
- Irrigation bottle (400mL)
- Foot control

RRP \$9,285

PROMO PRICE \$7,892

SAVE \$1,393



MODEL VARIOS COMBI PRO PERIO SET
ORDER CODE Y1003042

Additional to the Basic set, the Perio Set allows safe subgingival treatment for a range of procedures.

Includes:

- PERIO-CHAMBER
- PERIO-Hose
- PERIO-Handpiece
- Perio-Mate Nozzle tip
- Perio-Mate Nozzle tip remover
- Perio-Mate powder

RRP \$3,356

PROMO PRICE \$2,853

SAVE \$503



For more information scan the QR code.



ENDO-MATE DT2

ENDODONTIC MOTOR

20% OFF ENDO-MATE DT2

6 head angles for ease of treatment

Depending on the position of treatment site, the head angle can be adjusted to the best angle.

- Power supply: AC 230 V 50/60 Hz and Battery pack
- Number of Programs: 9
- Max Torque: 7 Ncm (when using 20:1 head)
- Speed range: 100-13,000 min⁻¹ (with 20:1, 1:1 geared heads)
- Charge/Continuous use: 5 h / 2 h (depending on workload)
- Dimensions: W 93 x D 147 x H 100 mm
- Weight: Unit 440 g/Motor handpiece and Cord 95 g.



MODEL ENDO-MATE DT2 (230V) (WITH HEAD)
ORDER CODE Y1004230

Includes:

- Control unit
- Motor and Cord
- MP-F20R head (20:1)
- Handpiece stand
- AC Adaptor

RRP \$2,775

PROMO PRICE \$2,220

SAVE \$555

For more information scan the QR code.



iCare HANDPIECE MAINTENANCE AND CARE

25% OFF iCARE AND ADAPTORS

Automatically lubricate

iCare sets a new standard for handpiece maintenance systems. In one stylish and compact body, it offers you both ease of use and multiple functions to achieve complete lubrication. The proper solution volume is preset based on handpiece type, and maintenance can be initiated at the touch of a button. The chuck cleaning system keeps chucks pristine, and the purging select mode lets you alter purge time depending on the situation for truly comprehensive maintenance. Thanks to its stylish and compact body, iCare is also space efficient and totally in harmony with your working environment.

iCare complete sets

- Power supply: AC 100-240 V 50/60 Hz
- Air pressure: 0.5-0.6 MPa (5.0-6.0 kgf/cm²)
- Weight: 6.5 kg (C2 type/C3 type)
- Dimension: W 280 x D 240 x H 350 (mm)
- Solution tank capacity: 1.2 L
- Connects only to Midwest 4-hole port
- Air piping is necessary (5.0-6.0 kgf/cm²)

Includes:

- iCare unit
- Power cord
- Chuck cleaning nozzle
- Maintenance oil for iCare (1 litre)
- Mist filter set
- Oil absorber sheet
- Air tube



Dual Standard Connectors + Dual Rotation Connectors

MODEL **iCARE C2 TYPE**
ORDER CODE **Y1002787**

RRP \$6,487

PROMO PRICE \$4,865

SAVE \$1,622



Single Standard Connector + Triple Rotation Connectors

MODEL **iCARE C3 TYPE**
ORDER CODE **Y1002788**

RRP \$6,929

PROMO PRICE \$5,197

SAVE \$1,732

Adaptors available to suit NSK, KaVo, Sirona, W&H and Morita couplings.

For more information scan the QR code.



VOLVERE i7 HIGH TORQUE, HIGH EFFICIENCY CORELESS MOTOR

20% OFF VOLVERE i7

Light and Compact to Fit into the Smallest Space

The VOLVERE i7 is another advance in NSK's quest to provide dental technicians access to its globally renowned dental laboratory engines. The stylish VOLVERE i7 delivers a compact, advanced design, all essential features and excellent value for money.

Standard Type Motor Handpiece

- Speed: 1,000-35,000 min⁻¹
- Max. Torque: 4.1 Ncm
- Dimensions: L 145 x ø24.5 (mm)
- Weight: 192 g (without cord)
- Cord Length: 1.2 m

E-Type Motor Handpiece

- Speed: 1,000-35,000 min⁻¹
- Max. Torque: 4.1 Ncm
- Dimensions: L 109.3 x ø24.5 (mm)
- Weight: 118 g (without cord)
- Cord Length: 1.2 m



Standard Micromotor



E-Type Micromotor

Complete with Standard Micromotor

MODEL **VOLVERE i7 RM**
ORDER CODE **Y1003493**

RRP \$2,500

PROMO PRICE \$2,000

SAVE \$500

E-Type Micromotor Without Handpiece

MODEL **VOLVERE i7 E**
ORDER CODE **Y1003494**

RRP \$2,058

PROMO PRICE \$1,646

SAVE \$412



For more information scan the QR code.



Choose from a wide range of configurations to suit your needs, options available for both Standard or E-Type Micromotors.

All prices as listed. Valid until 30 June 2025. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in NZD and exclude GST.

YOUR LOCAL NSK BUSINESS DEVELOPMENT AND EDUCATION SPECIALISTS

NAME	TITLE	REGION	PHONE NUMBER	EMAIL
Gina Narbey	National Sales Manager	New Zealand	027 224 5804	gina@nskoceania.com.au
Jennika Charlton	Business Development Manager	Upper New Zealand	027 263 8327	jennika@nskoceania.com.au
Neralee Montague	Business Development Manager	Lower New Zealand	027 263 8337	neralee@nskoceania.com.au